

NECTAR[®] MEDICAL




syntrax[®]
NECTAR[®] MEDICAL
 NATIVE WHEY PROTEIN ISOLATE


PROTEIN	SUGAR	FAT	SODIUM
10 GRAMS	0 GRAMS	0 GRAMS	25 MILLIGRAMS

To use the detachable, reusable clip included in every bag of Nectar[®], follow the instructions below:

1. Cut off the clip with scissors or snap off by pulling the clip straight up.



2. Fold bag over at least twice and slide clip on the top for a tight lock.



Pro Tip: Clip is reusable for other products such as cereal, snacks, etc.

NET WT:
2lb (907g)

UNFLAVORED

Mixes Instantly

Contains rBST-Free, Grass-Fed Whey

Kosher & Halal

100% Natural Ingredients

Lactose & Gluten Free



Every Nectar[®] contains a FreshShield™ sachet which removes all oxygen inside the package leaving only inert gases, such as nitrogen. Without oxygen, these inert gases envelop and protect the delicate protein powder and flavoring molecules from oxidation, rancidity, and bacterial growth. Our proteins are unique in their ability to have the same quality and taste from the day of manufacturing to the "best by" date.

NECTAR[®] MEDICAL

syntrax[®]
NECTAR[®] MEDICAL
 NATIVE WHEY PROTEIN ISOLATE

Nectar[®] Medical is a 100% natural protein formula designed to exceed even the most stringent and demanding requirements of the medical community. It is designed around the highest-quality native whey protein isolate that derives from rBST-free cows that feed on grass and forage. Whey protein isolate is considered by many experts to be the single highest-quality source of protein for human nutrition. Its bioavailability is without compare, and its amino acids are utilized more efficiently than any other source of protein.

Because Nectar[®] Medical has a bland taste and is instantly soluble and dispersible, it can be mixed unnoticeably and easily into nearly any beverage or food item, such as milk, tea, juice, flavored water, cereal, soup, pudding or mashed potatoes. Likewise, Nectar[®] Medical is an ideal product to be used in tube-feeding applications due to its ease of use and health-promoting benefits.

Nectar[®] Medical is specifically formulated to be used as a medical food for renal (dialysis) patients. Meeting the distinctive needs of these individuals, Nectar[®] Medical provides a superior protein source with low levels of potassium and phosphorus.

Nectar[®] Medical is specifically formulated to be used as a medical food for bariatric surgery patients. Meeting the distinctive needs of these individuals, Nectar[®] Medical provides an easily digestible, superior protein source with low levels of fat, sugar and total carbohydrates.

Nutrition Facts	
82 servings per container	
Serving size	1/2 Scoop (11g)
Amount Per Serving	
Calories	40
	% Daily Value*
Total Fat 0g	0%*
Saturated Fat 0g	0%*
Trans Fat 0g	
Cholesterol 0mg	0%*
Sodium 25mg	1%*
Total Carbohydrate 0g	0%*
Dietary Fiber 0g	0%*
Total Sugars 0g	
Includes 0g Added Sugars	0%*
Protein 10g	
Vitamin D 0mcg	0%
Calcium 50mg	4%
Iron 0mg	0%
Potassium 30mg	0%

* The % Daily Value (DV) tells you how much a nutrient in a serving of food contributes to a daily diet. 2,000 calories a day is used for general nutrition advice.

Important Amino Acids Per 100 Grams Of Protein	
Arginine+	2.4g
Glutamine/Glutamic Acid+	7.4g
Histidine+	1.7g
Isoleucine*	6.7g
Leucine*	11.9g
Lysine+	9.7g
Methionine+	2.0g
Phenylalanine+	3.3g
Threonine+	7.2g
Tryptophan+	2.1g
Valine*	6.2g

+ Essential Amino Acids
 * Branched-Chain Amino Acids
 • Important Nonessential Amino Acids

INGREDIENTS: Whey protein isolate*, lecithin.
CONTAINS: Milk, Soy. Produced in a facility that uses egg, wheat and peanut products.
 • Ultrafiltered and undenatured; includes beta lactoglobulin, alpha lactalbumin, glycomacropeptides, immunoglobulin, b.s. albumin, protease peptone, lactoferrin, lactoperoxidase.

DIRECTIONS FOR ORAL USE: Based on protein requirements, vigorously mix one to two scoops of Nectar[®] Medical per 4–16 oz of your favorite beverage or food until thoroughly blended. When adding to foods or warm liquids, add protein slowly while stirring. Do not add to liquids above 130°F to prevent denaturing of protein.
DIRECTIONS FOR TUBE-FEEDING USE: Mix one scoop per 2–6 oz of water until desired consistency is achieved. Stir until completely dissolved. Administer by syringe through feeding tube.
 Notice: Use this product as a food supplement only. Do not use for weight reduction.
 Manufactured by:
 SIGB, Inc.
 PO Box 1715
 Cape Girardeau, MO 63702 USA
 (866) 333-5103 (7483)
 www.syntrax.com
 Made in the USA
 facebook.com/OfficialSyntrax
 instagram.com/OfficialSyntrax

